

PRODUCT CATALOG

Pediatric dental crowns for primary teeth

What we make, how it is packaged, and how to keep your
operatory stocked.



Two materials, one sizing system

KTR makes crowns for restoring primary teeth. Zirconia covers every primary tooth, stainless steel covers the molars, and both are sized and labeled the same way.



Zirconia



Stainless steel

WHAT IT IS	1200 MPa zirconia, tooth-colored and metal-free	Medical-grade stainless steel, pre-trimmed and pre-crimped
TEETH COVERED	Incisors, cuspids and molars	All eight primary molars
SIZES	Six per tooth, 102 tooth and size combinations	Six per molar, 48 tooth and size combinations
FITS	Standard, Slim and Mini	One, with margins you can adapt
TO GET STARTED	Anterior, posterior or lower incisor kit	Starter kit of 48
TO RESTOCK	Single crowns	Boxes of 5, or all 48 boxes as a set

How a KTR manufacturer code reads

Each slot in a KTR case is labeled with the tooth, in Universal and Palmer notation, and the size. The manufacturer code for the crown that goes in it is built from the same parts.

KTR-Zirc-Slim-E3-AUR3

Slim fit zirconia crown, upper right central incisor, size 3

Zirc-Slim

Material and fit.
Stainless steel refill boxes start KTR-SSC-RB.

E3

Universal tooth letter and crown size

AUR3

The same tooth in Palmer notation (A, upper right) and the size again

Start with a kit

A kit puts one of every size in a labeled case, so the right crown is on hand at the first appointment and you only reorder what you use.



Zirconia anterior kit

Upper central and lateral incisors and all four cuspids, sizes 1 to 6.

48 white crowns + 48 colored try-on crowns

KTR-Zirc-SK-A-[FIT]



Zirconia posterior kit

All eight primary molars, sizes 2 to 7.

48 white crowns + 48 colored try-on crowns

KTR-Zirc-SK-P-[FIT]



Zirconia lower incisor kit

One universal crown shape for teeth N, O, P and Q, sizes 1 to 6.

6 white crowns + 6 colored try-on crowns

KTR-Zirc-SK-U-[FIT]



Stainless steel starter kit

All eight primary molars, sizes 2 to 7.

48 crowns

KTR-SSC-SK

Try-on crowns

Each kit pairs every white crown with a colored try-on crown of the same size, also zirconia. Check the size with the try-on, which is reusable and autoclave-resistant, and the crown you cement stays clean.

Pro-Kit


If you would rather have more white crowns than try-ons, order any zirconia kit as a Pro-Kit: the try-on crowns are replaced with a second white crown of each size. Add -Pro to the kit code.

Fit

Zirconia kits come in Standard, Slim or Mini fit, and the fit name ends the code. The next page shows how they differ.

Zirconia crowns, one at a time

When a slot runs low, order that crown on its own. Its code is the fit prefix below plus the code from the table. Every size also comes as a colored zirconia try-on crown.



Standard

Standard mesial-distal width and occlusal-cervical height.

KTR-Zirc-Standard-E3-AUR3

Slim

1 mm narrower mesial-distally, same height.

KTR-Zirc-Slim-E3-AUR3

Mini

1 mm shorter occlusal-cervically, same width.

KTR-Zirc-Mini-E3-AUR3

INCISORS

PRIMARY TOOTH	SIZE 1	SIZE 2	SIZE 3	SIZE 4	SIZE 5	SIZE 6
D Upper Right Lateral Incisor	D1-BUR1	D2-BUR2	D3-BUR3	D4-BUR4	D5-BUR5	D6-BUR6
E Upper Right Central Incisor	E1-AUR1	E2-AUR2	E3-AUR3	E4-AUR4	E5-AUR5	E6-AUR6
F Upper Left Central Incisor	F1-AUL1	F2-AUL2	F3-AUL3	F4-AUL4	F5-AUL5	F6-AUL6
G Upper Left Lateral Incisor	G1-BUL1	G2-BUL2	G3-BUL3	G4-BUL4	G5-BUL5	G6-BUL6
N-Q Universal Lower Incisors	U1	U2	U3	U4	U5	U6

CUSPIDS

PRIMARY TOOTH	SIZE 1	SIZE 2	SIZE 3	SIZE 4	SIZE 5	SIZE 6
C Upper Right Cuspid	C1-CUR1	C2-CUR2	C3-CUR3	C4-CUR4	C5-CUR5	C6-CUR6
H Upper Left Cuspid	H1-CUL1	H2-CUL2	H3-CUL3	H4-CUL4	H5-CUL5	H6-CUL6
M Lower Left Cuspid	M1-CLL1	M2-CLL2	M3-CLL3	M4-CLL4	M5-CLL5	M6-CLL6
R Lower Right Cuspid	R1-CLR1	R2-CLR2	R3-CLR3	R4-CLR4	R5-CLR5	R6-CLR6

MOLARS

PRIMARY TOOTH	SIZE 2	SIZE 3	SIZE 4	SIZE 5	SIZE 6	SIZE 7
A Upper Right 2nd Molar	A2-EUR2	A3-EUR3	A4-EUR4	A5-EUR5	A6-EUR6	A7-EUR7
B Upper Right 1st Molar	B2-DUR2	B3-DUR3	B4-DUR4	B5-DUR5	B6-DUR6	B7-DUR7
I Upper Left 1st Molar	I2-DUL2	I3-DUL3	I4-DUL4	I5-DUL5	I6-DUL6	I7-DUL7
J Upper Left 2nd Molar	J2-EUL2	J3-EUL3	J4-EUL4	J5-EUL5	J6-EUL6	J7-EUL7
K Lower Left 2nd Molar	K2-ELL2	K3-ELL3	K4-ELL4	K5-ELL5	K6-ELL6	K7-ELL7
L Lower Left 1st Molar	L2-DLL2	L3-DLL3	L4-DLL4	L5-DLL5	L6-DLL6	L7-DLL7
S Lower Right 1st Molar	S2-DLR2	S3-DLR3	S4-DLR4	S5-DLR5	S6-DLR6	S7-DLR7
T Lower Right 2nd Molar	T2-ELR2	T3-ELR3	T4-ELR4	T5-ELR5	T6-ELR6	T7-ELR7

Incisor and cuspid crowns have internal retention grooves designed to enhance bond retention.

Stainless steel crowns, by the box

Stainless steel crowns arrive pre-trimmed, belled and crimped, with the size laser-etched on the crown. Restock a single size with a box of 5, or the whole range in one order.



REFILL BOXES OF 5: KTR-SSC-RB- THEN THE CODE, E.G. KTR-SSC-RB-A4-EUR4

PRIMARY TOOTH	SIZE 2	SIZE 3	SIZE 4	SIZE 5	SIZE 6	SIZE 7
A Upper Right 2nd Molar	A2-EUR2	A3-EUR3	A4-EUR4	A5-EUR5	A6-EUR6	A7-EUR7
B Upper Right 1st Molar	B2-DUR2	B3-DUR3	B4-DUR4	B5-DUR5	B6-DUR6	B7-DUR7
I Upper Left 1st Molar	I2-DUL2	I3-DUL3	I4-DUL4	I5-DUL5	I6-DUL6	I7-DUL7
J Upper Left 2nd Molar	J2-EUL2	J3-EUL3	J4-EUL4	J5-EUL5	J6-EUL6	J7-EUL7
K Lower Left 2nd Molar	K2-ELL2	K3-ELL3	K4-ELL4	K5-ELL5	K6-ELL6	K7-ELL7
L Lower Left 1st Molar	L2-DLL2	L3-DLL3	L4-DLL4	L5-DLL5	L6-DLL6	L7-DLL7
S Lower Right 1st Molar	S2-DLR2	S3-DLR3	S4-DLR4	S5-DLR5	S6-DLR6	S7-DLR7
T Lower Right 2nd Molar	T2-ELR2	T3-ELR3	T4-ELR4	T5-ELR5	T6-ELR6	T7-ELR7



48 box refill set

One box of every size, 240 crowns, in a single order.

KTR-SSC-RB-48S



Crimping pliers

For adjusting a margin when a crown needs further contouring.

KTR-SSC-CP

A case for crowns you already have

The storage case from the kits is sold empty, for a second operatory or for crowns bought one at a time. Slots are labeled in Universal and Palmer notation, and the lid locks for transport.

Incisors and cuspids case	Teeth C, D, E, F, G, H, M and R, sizes 1 to 6	KTR-Case-A
Molars case	Teeth A, B, I, J, K, L, S and T, sizes 2 to 7	KTR-Case-P



How to order

ONLINE

ktrcrowns.com/shop

BY PHONE

1-800-506-5108

BY EMAIL

info@ktrcrowns.com

THROUGH YOUR MARKETPLACE

Search for the KTR
manufacturer code.

Orders placed with KTR before 3:00 PM Pacific, Monday to Friday, ship the same day. Instructions for use for both crown lines are at ktrcrowns.com.

Orthodontic supplies

KTR also supplies brackets, archwires, elastomerics, bonding and instruments. They have their own catalog at ktrcrowns.com/shop/ortho.

How to read a KTR manufacturer code

KTR

Brand

Zirc

Material: Zirc or SSC

Slim

Fit: Standard, Slim, Mini

E3

Universal tooth E, size 3

AUR3

Palmer A, upper right, size 3

ZIRCONIA STARTER KITS 48 white + 48 colored try-on crowns

KTR-Zirc-SK-A-Standard	Anterior kit, Standard fit
KTR-Zirc-SK-P-Standard	Posterior kit, Standard fit
KTR-Zirc-SK-A-Slim	Anterior kit, Slim fit
KTR-Zirc-SK-P-Slim	Posterior kit, Slim fit
KTR-Zirc-SK-A-Mini	Anterior kit, Mini fit
KTR-Zirc-SK-P-Mini	Posterior kit, Mini fit

ZIRCONIA CROWNS, STANDARD FIT 1 crown

KTR-Zirc-Standard-A2-EUR2	Upper Right 2nd Molar	Size 2
KTR-Zirc-Standard-A3-EUR3	Upper Right 2nd Molar	Size 3
KTR-Zirc-Standard-A4-EUR4	Upper Right 2nd Molar	Size 4
KTR-Zirc-Standard-A5-EUR5	Upper Right 2nd Molar	Size 5
KTR-Zirc-Standard-A6-EUR6	Upper Right 2nd Molar	Size 6
KTR-Zirc-Standard-A7-EUR7	Upper Right 2nd Molar	Size 7
KTR-Zirc-Standard-B2-DUR2	Upper Right 1st Molar	Size 2
KTR-Zirc-Standard-B3-DUR3	Upper Right 1st Molar	Size 3
KTR-Zirc-Standard-B4-DUR4	Upper Right 1st Molar	Size 4
KTR-Zirc-Standard-B5-DUR5	Upper Right 1st Molar	Size 5
KTR-Zirc-Standard-B6-DUR6	Upper Right 1st Molar	Size 6
KTR-Zirc-Standard-B7-DUR7	Upper Right 1st Molar	Size 7
KTR-Zirc-Standard-C1-CUR1	Upper Right Cuspid	Size 1
KTR-Zirc-Standard-C2-CUR2	Upper Right Cuspid	Size 2
KTR-Zirc-Standard-C3-CUR3	Upper Right Cuspid	Size 3
KTR-Zirc-Standard-C4-CUR4	Upper Right Cuspid	Size 4
KTR-Zirc-Standard-C5-CUR5	Upper Right Cuspid	Size 5
KTR-Zirc-Standard-C6-CUR6	Upper Right Cuspid	Size 6
KTR-Zirc-Standard-D1-BUR1	Upper Right Lateral Incisor	Size 1
KTR-Zirc-Standard-D2-BUR2	Upper Right Lateral Incisor	Size 2
KTR-Zirc-Standard-D3-BUR3	Upper Right Lateral Incisor	Size 3
KTR-Zirc-Standard-D4-BUR4	Upper Right Lateral Incisor	Size 4
KTR-Zirc-Standard-D5-BUR5	Upper Right Lateral Incisor	Size 5
KTR-Zirc-Standard-D6-BUR6	Upper Right Lateral Incisor	Size 6
KTR-Zirc-Standard-E1-AUR1	Upper Right Central Incisor	Size 1
KTR-Zirc-Standard-E2-AUR2	Upper Right Central Incisor	Size 2
KTR-Zirc-Standard-E3-AUR3	Upper Right Central Incisor	Size 3
KTR-Zirc-Standard-E4-AUR4	Upper Right Central Incisor	Size 4
KTR-Zirc-Standard-E5-AUR5	Upper Right Central Incisor	Size 5
KTR-Zirc-Standard-E6-AUR6	Upper Right Central Incisor	Size 6
KTR-Zirc-Standard-F1-AUL1	Upper Left Central Incisor	Size 1
KTR-Zirc-Standard-F2-AUL2	Upper Left Central Incisor	Size 2
KTR-Zirc-Standard-F3-AUL3	Upper Left Central Incisor	Size 3
KTR-Zirc-Standard-F4-AUL4	Upper Left Central Incisor	Size 4
KTR-Zirc-Standard-F5-AUL5	Upper Left Central Incisor	Size 5
KTR-Zirc-Standard-F6-AUL6	Upper Left Central Incisor	Size 6
KTR-Zirc-Standard-G1-BUL1	Upper Left Lateral Incisor	Size 1

ZIRCONIA CROWNS, STANDARD FIT continued

KTR-Zirc-Standard-G2-BUL2	Upper Left Lateral Incisor	Size 2
KTR-Zirc-Standard-G3-BUL3	Upper Left Lateral Incisor	Size 3
KTR-Zirc-Standard-G4-BUL4	Upper Left Lateral Incisor	Size 4
KTR-Zirc-Standard-G5-BUL5	Upper Left Lateral Incisor	Size 5
KTR-Zirc-Standard-G6-BUL6	Upper Left Lateral Incisor	Size 6
KTR-Zirc-Standard-H1-CUL1	Upper Left Cuspid	Size 1
KTR-Zirc-Standard-H2-CUL2	Upper Left Cuspid	Size 2
KTR-Zirc-Standard-H3-CUL3	Upper Left Cuspid	Size 3
KTR-Zirc-Standard-H4-CUL4	Upper Left Cuspid	Size 4
KTR-Zirc-Standard-H5-CUL5	Upper Left Cuspid	Size 5
KTR-Zirc-Standard-H6-CUL6	Upper Left Cuspid	Size 6
KTR-Zirc-Standard-I2-DUL2	Upper Left 1st Molar	Size 2
KTR-Zirc-Standard-I3-DUL3	Upper Left 1st Molar	Size 3
KTR-Zirc-Standard-I4-DUL4	Upper Left 1st Molar	Size 4
KTR-Zirc-Standard-I5-DUL5	Upper Left 1st Molar	Size 5
KTR-Zirc-Standard-I6-DUL6	Upper Left 1st Molar	Size 6
KTR-Zirc-Standard-I7-DUL7	Upper Left 1st Molar	Size 7
KTR-Zirc-Standard-J2-EUL2	Upper Left 2nd Molar	Size 2
KTR-Zirc-Standard-J3-EUL3	Upper Left 2nd Molar	Size 3
KTR-Zirc-Standard-J4-EUL4	Upper Left 2nd Molar	Size 4
KTR-Zirc-Standard-J5-EUL5	Upper Left 2nd Molar	Size 5
KTR-Zirc-Standard-J6-EUL6	Upper Left 2nd Molar	Size 6
KTR-Zirc-Standard-J7-EUL7	Upper Left 2nd Molar	Size 7
KTR-Zirc-Standard-K2-ELL2	Lower Left 2nd Molar	Size 2
KTR-Zirc-Standard-K3-ELL3	Lower Left 2nd Molar	Size 3
KTR-Zirc-Standard-K4-ELL4	Lower Left 2nd Molar	Size 4
KTR-Zirc-Standard-K5-ELL5	Lower Left 2nd Molar	Size 5
KTR-Zirc-Standard-K6-ELL6	Lower Left 2nd Molar	Size 6
KTR-Zirc-Standard-K7-ELL7	Lower Left 2nd Molar	Size 7
KTR-Zirc-Standard-L2-DLL2	Lower Left 1st Molar	Size 2
KTR-Zirc-Standard-L3-DLL3	Lower Left 1st Molar	Size 3
KTR-Zirc-Standard-L4-DLL4	Lower Left 1st Molar	Size 4
KTR-Zirc-Standard-L5-DLL5	Lower Left 1st Molar	Size 5
KTR-Zirc-Standard-L6-DLL6	Lower Left 1st Molar	Size 6
KTR-Zirc-Standard-L7-DLL7	Lower Left 1st Molar	Size 7
KTR-Zirc-Standard-M1-CLL1	Lower Left Cuspid	Size 1
KTR-Zirc-Standard-M2-CLL2	Lower Left Cuspid	Size 2
KTR-Zirc-Standard-M3-CLL3	Lower Left Cuspid	Size 3
KTR-Zirc-Standard-M4-CLL4	Lower Left Cuspid	Size 4
KTR-Zirc-Standard-M5-CLL5	Lower Left Cuspid	Size 5
KTR-Zirc-Standard-M6-CLL6	Lower Left Cuspid	Size 6
KTR-Zirc-Standard-R1-CLR1	Lower Right Cuspid	Size 1
KTR-Zirc-Standard-R2-CLR2	Lower Right Cuspid	Size 2
KTR-Zirc-Standard-R3-CLR3	Lower Right Cuspid	Size 3
KTR-Zirc-Standard-R4-CLR4	Lower Right Cuspid	Size 4

ZIRCONIA CROWNS, STANDARD FIT continued

KTR-Zirc-Standard-R5-CLR5	Lower Right Cuspid	Size 5
KTR-Zirc-Standard-R6-CLR6	Lower Right Cuspid	Size 6
KTR-Zirc-Standard-S2-DLR2	Lower Right 1st Molar	Size 2
KTR-Zirc-Standard-S3-DLR3	Lower Right 1st Molar	Size 3
KTR-Zirc-Standard-S4-DLR4	Lower Right 1st Molar	Size 4
KTR-Zirc-Standard-S5-DLR5	Lower Right 1st Molar	Size 5
KTR-Zirc-Standard-S6-DLR6	Lower Right 1st Molar	Size 6
KTR-Zirc-Standard-S7-DLR7	Lower Right 1st Molar	Size 7
KTR-Zirc-Standard-T2-ELR2	Lower Right 2nd Molar	Size 2
KTR-Zirc-Standard-T3-ELR3	Lower Right 2nd Molar	Size 3
KTR-Zirc-Standard-T4-ELR4	Lower Right 2nd Molar	Size 4
KTR-Zirc-Standard-T5-ELR5	Lower Right 2nd Molar	Size 5
KTR-Zirc-Standard-T6-ELR6	Lower Right 2nd Molar	Size 6
KTR-Zirc-Standard-T7-ELR7	Lower Right 2nd Molar	Size 7
KTR-Zirc-Standard-U1	Universal Lower Incisors	Size 1
KTR-Zirc-Standard-U2	Universal Lower Incisors	Size 2
KTR-Zirc-Standard-U3	Universal Lower Incisors	Size 3
KTR-Zirc-Standard-U4	Universal Lower Incisors	Size 4
KTR-Zirc-Standard-U5	Universal Lower Incisors	Size 5
KTR-Zirc-Standard-U6	Universal Lower Incisors	Size 6

ZIRCONIA CROWNS, SLIM FIT 1 crown

KTR-Zirc-Slim-A2-EUR2	Upper Right 2nd Molar	Size 2
KTR-Zirc-Slim-A3-EUR3	Upper Right 2nd Molar	Size 3
KTR-Zirc-Slim-A4-EUR4	Upper Right 2nd Molar	Size 4
KTR-Zirc-Slim-A5-EUR5	Upper Right 2nd Molar	Size 5
KTR-Zirc-Slim-A6-EUR6	Upper Right 2nd Molar	Size 6
KTR-Zirc-Slim-A7-EUR7	Upper Right 2nd Molar	Size 7
KTR-Zirc-Slim-B2-DUR2	Upper Right 1st Molar	Size 2
KTR-Zirc-Slim-B3-DUR3	Upper Right 1st Molar	Size 3
KTR-Zirc-Slim-B4-DUR4	Upper Right 1st Molar	Size 4
KTR-Zirc-Slim-B5-DUR5	Upper Right 1st Molar	Size 5
KTR-Zirc-Slim-B6-DUR6	Upper Right 1st Molar	Size 6
KTR-Zirc-Slim-B7-DUR7	Upper Right 1st Molar	Size 7
KTR-Zirc-Slim-C1-CUR1	Upper Right Cuspid	Size 1
KTR-Zirc-Slim-C2-CUR2	Upper Right Cuspid	Size 2
KTR-Zirc-Slim-C3-CUR3	Upper Right Cuspid	Size 3
KTR-Zirc-Slim-C4-CUR4	Upper Right Cuspid	Size 4
KTR-Zirc-Slim-C5-CUR5	Upper Right Cuspid	Size 5
KTR-Zirc-Slim-C6-CUR6	Upper Right Cuspid	Size 6
KTR-Zirc-Slim-D1-BUR1	Upper Right Lateral Incisor	Size 1
KTR-Zirc-Slim-D2-BUR2	Upper Right Lateral Incisor	Size 2
KTR-Zirc-Slim-D3-BUR3	Upper Right Lateral Incisor	Size 3
KTR-Zirc-Slim-D4-BUR4	Upper Right Lateral Incisor	Size 4
KTR-Zirc-Slim-D5-BUR5	Upper Right Lateral Incisor	Size 5
KTR-Zirc-Slim-D6-BUR6	Upper Right Lateral Incisor	Size 6
KTR-Zirc-Slim-E1-AUR1	Upper Right Central Incisor	Size 1
KTR-Zirc-Slim-E2-AUR2	Upper Right Central Incisor	Size 2
KTR-Zirc-Slim-E3-AUR3	Upper Right Central Incisor	Size 3
KTR-Zirc-Slim-E4-AUR4	Upper Right Central Incisor	Size 4

ZIRCONIA CROWNS, SLIM FIT continued

KTR-Zirc-Slim-E5-AUR5	Upper Right Central Incisor	Size 5
KTR-Zirc-Slim-E6-AUR6	Upper Right Central Incisor	Size 6
KTR-Zirc-Slim-F1-AUL1	Upper Left Central Incisor	Size 1
KTR-Zirc-Slim-F2-AUL2	Upper Left Central Incisor	Size 2
KTR-Zirc-Slim-F3-AUL3	Upper Left Central Incisor	Size 3
KTR-Zirc-Slim-F4-AUL4	Upper Left Central Incisor	Size 4
KTR-Zirc-Slim-F5-AUL5	Upper Left Central Incisor	Size 5
KTR-Zirc-Slim-F6-AUL6	Upper Left Central Incisor	Size 6
KTR-Zirc-Slim-G1-BUL1	Upper Left Lateral Incisor	Size 1
KTR-Zirc-Slim-G2-BUL2	Upper Left Lateral Incisor	Size 2
KTR-Zirc-Slim-G3-BUL3	Upper Left Lateral Incisor	Size 3
KTR-Zirc-Slim-G4-BUL4	Upper Left Lateral Incisor	Size 4
KTR-Zirc-Slim-G5-BUL5	Upper Left Lateral Incisor	Size 5
KTR-Zirc-Slim-G6-BUL6	Upper Left Lateral Incisor	Size 6
KTR-Zirc-Slim-H1-CUL1	Upper Left Cuspid	Size 1
KTR-Zirc-Slim-H2-CUL2	Upper Left Cuspid	Size 2
KTR-Zirc-Slim-H3-CUL3	Upper Left Cuspid	Size 3
KTR-Zirc-Slim-H4-CUL4	Upper Left Cuspid	Size 4
KTR-Zirc-Slim-H5-CUL5	Upper Left Cuspid	Size 5
KTR-Zirc-Slim-H6-CUL6	Upper Left Cuspid	Size 6
KTR-Zirc-Slim-I2-DUL2	Upper Left 1st Molar	Size 2
KTR-Zirc-Slim-I3-DUL3	Upper Left 1st Molar	Size 3
KTR-Zirc-Slim-I4-DUL4	Upper Left 1st Molar	Size 4
KTR-Zirc-Slim-I5-DUL5	Upper Left 1st Molar	Size 5
KTR-Zirc-Slim-I6-DUL6	Upper Left 1st Molar	Size 6
KTR-Zirc-Slim-I7-DUL7	Upper Left 1st Molar	Size 7
KTR-Zirc-Slim-J2-EUL2	Upper Left 2nd Molar	Size 2
KTR-Zirc-Slim-J3-EUL3	Upper Left 2nd Molar	Size 3
KTR-Zirc-Slim-J4-EUL4	Upper Left 2nd Molar	Size 4
KTR-Zirc-Slim-J5-EUL5	Upper Left 2nd Molar	Size 5
KTR-Zirc-Slim-J6-EUL6	Upper Left 2nd Molar	Size 6
KTR-Zirc-Slim-J7-EUL7	Upper Left 2nd Molar	Size 7
KTR-Zirc-Slim-K2-ELL2	Lower Left 2nd Molar	Size 2
KTR-Zirc-Slim-K3-ELL3	Lower Left 2nd Molar	Size 3
KTR-Zirc-Slim-K4-ELL4	Lower Left 2nd Molar	Size 4
KTR-Zirc-Slim-K5-ELL5	Lower Left 2nd Molar	Size 5
KTR-Zirc-Slim-K6-ELL6	Lower Left 2nd Molar	Size 6
KTR-Zirc-Slim-K7-ELL7	Lower Left 2nd Molar	Size 7
KTR-Zirc-Slim-L2-DLL2	Lower Left 1st Molar	Size 2
KTR-Zirc-Slim-L3-DLL3	Lower Left 1st Molar	Size 3
KTR-Zirc-Slim-L4-DLL4	Lower Left 1st Molar	Size 4
KTR-Zirc-Slim-L5-DLL5	Lower Left 1st Molar	Size 5
KTR-Zirc-Slim-L6-DLL6	Lower Left 1st Molar	Size 6
KTR-Zirc-Slim-L7-DLL7	Lower Left 1st Molar	Size 7
KTR-Zirc-Slim-M1-CLL1	Lower Left Cuspid	Size 1
KTR-Zirc-Slim-M2-CLL2	Lower Left Cuspid	Size 2
KTR-Zirc-Slim-M3-CLL3	Lower Left Cuspid	Size 3
KTR-Zirc-Slim-M4-CLL4	Lower Left Cuspid	Size 4
KTR-Zirc-Slim-M5-CLL5	Lower Left Cuspid	Size 5
KTR-Zirc-Slim-M6-CLL6	Lower Left Cuspid	Size 6

ZIRCONIA CROWNS, SLIM FIT continued

KTR-Zirc-Slim-R1-CLR1	Lower Right Cuspid	Size 1
KTR-Zirc-Slim-R2-CLR2	Lower Right Cuspid	Size 2
KTR-Zirc-Slim-R3-CLR3	Lower Right Cuspid	Size 3
KTR-Zirc-Slim-R4-CLR4	Lower Right Cuspid	Size 4
KTR-Zirc-Slim-R5-CLR5	Lower Right Cuspid	Size 5
KTR-Zirc-Slim-R6-CLR6	Lower Right Cuspid	Size 6
KTR-Zirc-Slim-S2-DLR2	Lower Right 1st Molar	Size 2
KTR-Zirc-Slim-S3-DLR3	Lower Right 1st Molar	Size 3
KTR-Zirc-Slim-S4-DLR4	Lower Right 1st Molar	Size 4
KTR-Zirc-Slim-S5-DLR5	Lower Right 1st Molar	Size 5
KTR-Zirc-Slim-S6-DLR6	Lower Right 1st Molar	Size 6
KTR-Zirc-Slim-S7-DLR7	Lower Right 1st Molar	Size 7
KTR-Zirc-Slim-T2-ELR2	Lower Right 2nd Molar	Size 2
KTR-Zirc-Slim-T3-ELR3	Lower Right 2nd Molar	Size 3
KTR-Zirc-Slim-T4-ELR4	Lower Right 2nd Molar	Size 4
KTR-Zirc-Slim-T5-ELR5	Lower Right 2nd Molar	Size 5
KTR-Zirc-Slim-T6-ELR6	Lower Right 2nd Molar	Size 6
KTR-Zirc-Slim-T7-ELR7	Lower Right 2nd Molar	Size 7
KTR-Zirc-Slim-U1	Universal Lower Incisors	Size 1
KTR-Zirc-Slim-U2	Universal Lower Incisors	Size 2
KTR-Zirc-Slim-U3	Universal Lower Incisors	Size 3
KTR-Zirc-Slim-U4	Universal Lower Incisors	Size 4
KTR-Zirc-Slim-U5	Universal Lower Incisors	Size 5
KTR-Zirc-Slim-U6	Universal Lower Incisors	Size 6

ZIRCONIA CROWNS, MINI FIT 1 crown

KTR-Zirc-Mini-A2-EUR2	Upper Right 2nd Molar	Size 2
KTR-Zirc-Mini-A3-EUR3	Upper Right 2nd Molar	Size 3
KTR-Zirc-Mini-A4-EUR4	Upper Right 2nd Molar	Size 4
KTR-Zirc-Mini-A5-EUR5	Upper Right 2nd Molar	Size 5
KTR-Zirc-Mini-A6-EUR6	Upper Right 2nd Molar	Size 6
KTR-Zirc-Mini-A7-EUR7	Upper Right 2nd Molar	Size 7
KTR-Zirc-Mini-B2-DUR2	Upper Right 1st Molar	Size 2
KTR-Zirc-Mini-B3-DUR3	Upper Right 1st Molar	Size 3
KTR-Zirc-Mini-B4-DUR4	Upper Right 1st Molar	Size 4
KTR-Zirc-Mini-B5-DUR5	Upper Right 1st Molar	Size 5
KTR-Zirc-Mini-B6-DUR6	Upper Right 1st Molar	Size 6
KTR-Zirc-Mini-B7-DUR7	Upper Right 1st Molar	Size 7
KTR-Zirc-Mini-C1-CUR1	Upper Right Cuspid	Size 1
KTR-Zirc-Mini-C2-CUR2	Upper Right Cuspid	Size 2
KTR-Zirc-Mini-C3-CUR3	Upper Right Cuspid	Size 3
KTR-Zirc-Mini-C4-CUR4	Upper Right Cuspid	Size 4
KTR-Zirc-Mini-C5-CUR5	Upper Right Cuspid	Size 5
KTR-Zirc-Mini-C6-CUR6	Upper Right Cuspid	Size 6
KTR-Zirc-Mini-D1-BUR1	Upper Right Lateral Incisor	Size 1
KTR-Zirc-Mini-D2-BUR2	Upper Right Lateral Incisor	Size 2
KTR-Zirc-Mini-D3-BUR3	Upper Right Lateral Incisor	Size 3
KTR-Zirc-Mini-D4-BUR4	Upper Right Lateral Incisor	Size 4
KTR-Zirc-Mini-D5-BUR5	Upper Right Lateral Incisor	Size 5
KTR-Zirc-Mini-D6-BUR6	Upper Right Lateral Incisor	Size 6

ZIRCONIA CROWNS, MINI FIT continued

KTR-Zirc-Mini-E1-AUR1	Upper Right Central Incisor	Size 1
KTR-Zirc-Mini-E2-AUR2	Upper Right Central Incisor	Size 2
KTR-Zirc-Mini-E3-AUR3	Upper Right Central Incisor	Size 3
KTR-Zirc-Mini-E4-AUR4	Upper Right Central Incisor	Size 4
KTR-Zirc-Mini-E5-AUR5	Upper Right Central Incisor	Size 5
KTR-Zirc-Mini-E6-AUR6	Upper Right Central Incisor	Size 6
KTR-Zirc-Mini-F1-AUL1	Upper Left Central Incisor	Size 1
KTR-Zirc-Mini-F2-AUL2	Upper Left Central Incisor	Size 2
KTR-Zirc-Mini-F3-AUL3	Upper Left Central Incisor	Size 3
KTR-Zirc-Mini-F4-AUL4	Upper Left Central Incisor	Size 4
KTR-Zirc-Mini-F5-AUL5	Upper Left Central Incisor	Size 5
KTR-Zirc-Mini-F6-AUL6	Upper Left Central Incisor	Size 6
KTR-Zirc-Mini-G1-BUL1	Upper Left Lateral Incisor	Size 1
KTR-Zirc-Mini-G2-BUL2	Upper Left Lateral Incisor	Size 2
KTR-Zirc-Mini-G3-BUL3	Upper Left Lateral Incisor	Size 3
KTR-Zirc-Mini-G4-BUL4	Upper Left Lateral Incisor	Size 4
KTR-Zirc-Mini-G5-BUL5	Upper Left Lateral Incisor	Size 5
KTR-Zirc-Mini-G6-BUL6	Upper Left Lateral Incisor	Size 6
KTR-Zirc-Mini-H1-CUL1	Upper Left Cuspid	Size 1
KTR-Zirc-Mini-H2-CUL2	Upper Left Cuspid	Size 2
KTR-Zirc-Mini-H3-CUL3	Upper Left Cuspid	Size 3
KTR-Zirc-Mini-H4-CUL4	Upper Left Cuspid	Size 4
KTR-Zirc-Mini-H5-CUL5	Upper Left Cuspid	Size 5
KTR-Zirc-Mini-H6-CUL6	Upper Left Cuspid	Size 6
KTR-Zirc-Mini-I2-DUL2	Upper Left 1st Molar	Size 2
KTR-Zirc-Mini-I3-DUL3	Upper Left 1st Molar	Size 3
KTR-Zirc-Mini-I4-DUL4	Upper Left 1st Molar	Size 4
KTR-Zirc-Mini-I5-DUL5	Upper Left 1st Molar	Size 5
KTR-Zirc-Mini-I6-DUL6	Upper Left 1st Molar	Size 6
KTR-Zirc-Mini-I7-DUL7	Upper Left 1st Molar	Size 7
KTR-Zirc-Mini-J2-EUL2	Upper Left 2nd Molar	Size 2
KTR-Zirc-Mini-J3-EUL3	Upper Left 2nd Molar	Size 3
KTR-Zirc-Mini-J4-EUL4	Upper Left 2nd Molar	Size 4
KTR-Zirc-Mini-J5-EUL5	Upper Left 2nd Molar	Size 5
KTR-Zirc-Mini-J6-EUL6	Upper Left 2nd Molar	Size 6
KTR-Zirc-Mini-J7-EUL7	Upper Left 2nd Molar	Size 7
KTR-Zirc-Mini-K2-ELL2	Lower Left 2nd Molar	Size 2
KTR-Zirc-Mini-K3-ELL3	Lower Left 2nd Molar	Size 3
KTR-Zirc-Mini-K4-ELL4	Lower Left 2nd Molar	Size 4
KTR-Zirc-Mini-K5-ELL5	Lower Left 2nd Molar	Size 5
KTR-Zirc-Mini-K6-ELL6	Lower Left 2nd Molar	Size 6
KTR-Zirc-Mini-K7-ELL7	Lower Left 2nd Molar	Size 7
KTR-Zirc-Mini-L2-DLL2	Lower Left 1st Molar	Size 2
KTR-Zirc-Mini-L3-DLL3	Lower Left 1st Molar	Size 3
KTR-Zirc-Mini-L4-DLL4	Lower Left 1st Molar	Size 4
KTR-Zirc-Mini-L5-DLL5	Lower Left 1st Molar	Size 5
KTR-Zirc-Mini-L6-DLL6	Lower Left 1st Molar	Size 6
KTR-Zirc-Mini-L7-DLL7	Lower Left 1st Molar	Size 7
KTR-Zirc-Mini-M1-CLL1	Lower Left Cuspid	Size 1
KTR-Zirc-Mini-M2-CLL2	Lower Left Cuspid	Size 2

ZIRCONIA CROWNS, MINI FIT continued

KTR-Zirc-Mini-M3-CLL3	Lower Left Cuspid	Size 3
KTR-Zirc-Mini-M4-CLL4	Lower Left Cuspid	Size 4
KTR-Zirc-Mini-M5-CLL5	Lower Left Cuspid	Size 5
KTR-Zirc-Mini-M6-CLL6	Lower Left Cuspid	Size 6
KTR-Zirc-Mini-R1-CLR1	Lower Right Cuspid	Size 1
KTR-Zirc-Mini-R2-CLR2	Lower Right Cuspid	Size 2
KTR-Zirc-Mini-R3-CLR3	Lower Right Cuspid	Size 3
KTR-Zirc-Mini-R4-CLR4	Lower Right Cuspid	Size 4
KTR-Zirc-Mini-R5-CLR5	Lower Right Cuspid	Size 5
KTR-Zirc-Mini-R6-CLR6	Lower Right Cuspid	Size 6
KTR-Zirc-Mini-S2-DLR2	Lower Right 1st Molar	Size 2
KTR-Zirc-Mini-S3-DLR3	Lower Right 1st Molar	Size 3
KTR-Zirc-Mini-S4-DLR4	Lower Right 1st Molar	Size 4
KTR-Zirc-Mini-S5-DLR5	Lower Right 1st Molar	Size 5
KTR-Zirc-Mini-S6-DLR6	Lower Right 1st Molar	Size 6
KTR-Zirc-Mini-S7-DLR7	Lower Right 1st Molar	Size 7
KTR-Zirc-Mini-T2-ELR2	Lower Right 2nd Molar	Size 2
KTR-Zirc-Mini-T3-ELR3	Lower Right 2nd Molar	Size 3
KTR-Zirc-Mini-T4-ELR4	Lower Right 2nd Molar	Size 4
KTR-Zirc-Mini-T5-ELR5	Lower Right 2nd Molar	Size 5
KTR-Zirc-Mini-T6-ELR6	Lower Right 2nd Molar	Size 6
KTR-Zirc-Mini-T7-ELR7	Lower Right 2nd Molar	Size 7
KTR-Zirc-Mini-U1	Universal Lower Incisors	Size 1
KTR-Zirc-Mini-U2	Universal Lower Incisors	Size 2
KTR-Zirc-Mini-U3	Universal Lower Incisors	Size 3
KTR-Zirc-Mini-U4	Universal Lower Incisors	Size 4
KTR-Zirc-Mini-U5	Universal Lower Incisors	Size 5
KTR-Zirc-Mini-U6	Universal Lower Incisors	Size 6

STAINLESS STEEL STARTER KIT 48 crowns

KTR-SSC-SK Starter kit, primary molars

STAINLESS STEEL REFILL BOXES 5 crowns per box

KTR-SSC-RB-A2-EUR2	Upper Right 2nd Molar	Size 2
KTR-SSC-RB-A3-EUR3	Upper Right 2nd Molar	Size 3
KTR-SSC-RB-A4-EUR4	Upper Right 2nd Molar	Size 4
KTR-SSC-RB-A5-EUR5	Upper Right 2nd Molar	Size 5
KTR-SSC-RB-A6-EUR6	Upper Right 2nd Molar	Size 6
KTR-SSC-RB-A7-EUR7	Upper Right 2nd Molar	Size 7
KTR-SSC-RB-B2-DUR2	Upper Right 1st Molar	Size 2
KTR-SSC-RB-B3-DUR3	Upper Right 1st Molar	Size 3
KTR-SSC-RB-B4-DUR4	Upper Right 1st Molar	Size 4
KTR-SSC-RB-B5-DUR5	Upper Right 1st Molar	Size 5
KTR-SSC-RB-B6-DUR6	Upper Right 1st Molar	Size 6
KTR-SSC-RB-B7-DUR7	Upper Right 1st Molar	Size 7
KTR-SSC-RB-I2-DUL2	Upper Left 1st Molar	Size 2
KTR-SSC-RB-I3-DUL3	Upper Left 1st Molar	Size 3
KTR-SSC-RB-I4-DUL4	Upper Left 1st Molar	Size 4
KTR-SSC-RB-I5-DUL5	Upper Left 1st Molar	Size 5
KTR-SSC-RB-I6-DUL6	Upper Left 1st Molar	Size 6
KTR-SSC-RB-I7-DUL7	Upper Left 1st Molar	Size 7

STAINLESS STEEL REFILL BOXES continued

KTR-SSC-RB-J2-EUL2	Upper Left 2nd Molar	Size 2
KTR-SSC-RB-J3-EUL3	Upper Left 2nd Molar	Size 3
KTR-SSC-RB-J4-EUL4	Upper Left 2nd Molar	Size 4
KTR-SSC-RB-J5-EUL5	Upper Left 2nd Molar	Size 5
KTR-SSC-RB-J6-EUL6	Upper Left 2nd Molar	Size 6
KTR-SSC-RB-J7-EUL7	Upper Left 2nd Molar	Size 7
KTR-SSC-RB-K2-ELL2	Lower Left 2nd Molar	Size 2
KTR-SSC-RB-K3-ELL3	Lower Left 2nd Molar	Size 3
KTR-SSC-RB-K4-ELL4	Lower Left 2nd Molar	Size 4
KTR-SSC-RB-K5-ELL5	Lower Left 2nd Molar	Size 5
KTR-SSC-RB-K6-ELL6	Lower Left 2nd Molar	Size 6
KTR-SSC-RB-K7-ELL7	Lower Left 2nd Molar	Size 7
KTR-SSC-RB-L2-DLL2	Lower Left 1st Molar	Size 2
KTR-SSC-RB-L3-DLL3	Lower Left 1st Molar	Size 3
KTR-SSC-RB-L4-DLL4	Lower Left 1st Molar	Size 4
KTR-SSC-RB-L5-DLL5	Lower Left 1st Molar	Size 5
KTR-SSC-RB-L6-DLL6	Lower Left 1st Molar	Size 6
KTR-SSC-RB-L7-DLL7	Lower Left 1st Molar	Size 7
KTR-SSC-RB-S2-DLR2	Lower Right 1st Molar	Size 2
KTR-SSC-RB-S3-DLR3	Lower Right 1st Molar	Size 3
KTR-SSC-RB-S4-DLR4	Lower Right 1st Molar	Size 4
KTR-SSC-RB-S5-DLR5	Lower Right 1st Molar	Size 5
KTR-SSC-RB-S6-DLR6	Lower Right 1st Molar	Size 6
KTR-SSC-RB-S7-DLR7	Lower Right 1st Molar	Size 7
KTR-SSC-RB-T2-ELR2	Lower Right 2nd Molar	Size 2
KTR-SSC-RB-T3-ELR3	Lower Right 2nd Molar	Size 3
KTR-SSC-RB-T4-ELR4	Lower Right 2nd Molar	Size 4
KTR-SSC-RB-T5-ELR5	Lower Right 2nd Molar	Size 5
KTR-SSC-RB-T6-ELR6	Lower Right 2nd Molar	Size 6
KTR-SSC-RB-T7-ELR7	Lower Right 2nd Molar	Size 7

STAINLESS STEEL REFILL SET 48 boxes, 240 crowns

KTR-SSC-RB-48S One refill box of every size

ACCESSORIES

KTR-Case-A	Storage case, incisors and cuspids
KTR-Case-P	Storage case, molars
KTR-SSC-CP	Crimping pliers

ORDER BY CODE ktrcrowns.com/shop 1-800-506-5108 info@ktrcrowns.com or search the code on your marketplace